

Intensify your skincare easy, effective & fast

IB Skin-Roller (Gold Plated) COSMETIC

patented, gold plated, durable, elegant, easy to clean

Short summary:

Handy, micro-needling gold roller with 192 needles to intensify the absorption of serums and gels.

Description:

Handy, micro-needling gold roller for non-surgical and non-ablative beauty treatments, mainly designed for the face and neck area. However, the roller can also be used for other parts of the body such as hands, feet, chest, arms and other areas.

Cosmetic Needling is a simple but effective technique to rejuvenate the upper layer of the skin. As there are neither nerves nor blood vessels in the epidermis, the treatment is painless. Clinical studies have shown that the absorption of serums or gels increases up to 100 times, when micro-needling is applied beforehand.



IB Skin-Rollers (Cosmetic) are suitable for personal use:

- Stimulates the skin
- Intensifies absorption of serums or gels
- Quickly regenerates the cells with appropriate skincare products
- Optimizes daily skincare products
- Promotes the absorption of nutrients into the skin
- Gold, ideal for sensitive skin

Product Information:

Material: brass | Head: 14 carat gold plated | Size: 13 cm
Weight: 107 g (including cap)

The IB Skin-Roller is designed to be used by one person only; for reasons of hygiene **it should not be shared**. Thoroughly clean the roller before and after each use. The 192 needles of the roller head should never be placed on hard surfaces. They must be handled with care and should not be damaged in any way. Pay attention to the contra-indications and recommendations. Keep the roller out of the reach of children. Use only as directed.

Application intensity:

After the skin is accustomed to the treatment the procedure can be carried out daily, preferably in the evening. In the beginning it is advisable to use the roller 1-2 times a week, before increasing frequency to daily use.

Storage: Store the appliance in a dry place.

Clean and dry the roller properly

Please clean the roller before and after each use. Rinse the roller under warm running water. Clean the roller once a week with a non-corrosive disinfectant. Do not put the roller away while it is still wet. Place the roller in the groove of the cap and let it dry completely, then place the protection cap over the roller.

Do not use in sensitive areas

When using the roller on sensitive areas such as the larynx, roll very carefully.

Apply appropriate skincare products only after micro-needling

After micro-needling, make sure the products you use are free from any fragrances and additives such as emulsifiers, mineral oils, colorants, preservatives or dyes, which can cause allergic reactions.