

Intensify your skincare easy, effective & fast

IB Skin-Stamp (Gold Plated) COSMETIC

patented, gold plated, durable, elegant, easy to clean

Short summary:

Handy, micro-needling gold stamp with 142 needles to intensify the absorption of serums and gels.

Description:

Handy, micro-needling gold stamp for non-surgical and non-ablative beauty treatments, mainly designed for the face and neck area. However, the stamp can also be used for other parts of the body such as hands, feet, chest, arms and other areas.

Cosmetic Needling is a simple but effective technique to rejuvenate the upper layer of the skin. As there are neither nerves nor blood vessels in the epidermis, the treatment is painless. Clinical studies have shown that the absorption of serums or gels increases up to 100 times, when micro-needling is applied beforehand.

IB Skin-Stamps (Cosmetic) are suitable for personal use:

- Stimulates the skin
- Intensifies absorption of serums or gels
- Quickly regenerates the cells with appropriate skincare products
- Optimizes daily skincare products
- Promotes the absorption of nutrients into the skin
- Gold, ideal for sensitive skin

Product Information:

Material: brass | Head: 14 carat gold plated | Size: 12 cm

Weight: 65 g (incl. cap)

The IB Skin-Stamp is designed to be used by one person only; for reasons of hygiene **it should not be shared**. Thoroughly clean the stamp before and after each use. The 142 needles of the stamp head should never be placed on hard surfaces. They must be handled with care and should not be damaged in any way. Pay attention to the contra-indications and recommendations. Keep the stamp out of the reach of children. Use only as directed.

APPLICATION INTENSITY (Cosmetic)

After the skin is accustomed to the treatment the procedure can be carried out daily, preferably in the evening. In the beginning it is advisable to use the stamp 1-2 times a week, before increasing the frequency to daily use.

Storage: Store the appliance in a dry place.

Clean and dry the stamp properly

Please clean the stamp before and after each use. Rinse the stamp under warm running water. Clean the stamp once a week with a non-corrosive disinfectant. Do not put the stamp away while it is still wet. Place the stamp in the groove of the cap and let it dry completely, then place the protection cap over the stamp.

Do not use in sensitive areas

When using the stamp on sensitive areas such as the larynx, be very carefully.

Apply appropriate skincare products only after micro-needling

After micro-needling, make sure the products you use are free from any fragrances and additives such as emulsifiers, mineral oils, colorants, preservatives or dyes, which can cause allergic reactions.

